



## Branran Medical Corporation

---

10015 Muirlands Road, Suite C • Irvine, CA 92618 • Tel: (949) 598-7166 • Fax: (949) 598-7167

### **Drug and Alcohol Testing Now in a Single, On-site Oral Fluid Device**

September 1, 2006

Branran Medical Corporation today announced a breakthrough in on-site drug and alcohol screening technology. The new product called OratectPlus™ combines drug and alcohol oral fluid testing in a single device. Never before has there been a test that can detect drugs plus alcohol in the same device. This new on-site saliva screen test allows test administrators to obtain the rapid results without the concern of handling a urine specimen.

By simply having the subject swab his/her mouth with the test device, drug and alcohol initial test results can be available in minutes. The reliability of an oral fluid test is comparable to urinalysis testing for the most commonly abused drugs.

OratectPlus™ simultaneously tests for 6 drugs including marijuana, cocaine, opiates, amphetamine, methamphetamine/ecstasy and either PCP or benzodiazepines plus alcohol. The test device consists of a collection pad on one end, two test panels in the middle for the drug results and a separate pad for the alcohol result. The test protocol entails swabbing the collection pad inside the mouth until blue lines in the test panels start to flow, indicating a sufficient amount of saliva has been collected. The presence of a red line corresponding to a specific drug on the panel indicates a negative result, while the absence of a red line indicates a presumptive positive result. A color change from white to green, blue or gray indicates a presumptive positive alcohol result. With a presumptive positive result, the collection pad is detached from the device and sent to a SAMHSA laboratory for further testing, allowing confirmation of the on-site result without the drawing of another sample.

Now Human Resource Managers and Criminal Justice Professionals who have the need for rapid, on-site drug and alcohol testing results but who have been reluctant to handle urine specimens have another option. OratectPlus™ may be the answer they have been looking for.

For more information on Branran Medical Corporation or its drug testing products, please visit [www.branranmedical.com](http://www.branranmedical.com).

#### About Branran Medical Corporation

Branran Medical Corporation, located in Irvine, California, is a manufacturer and marketer of rapid, on-site tests for drugs of abuse. Its products include adulteration tests, CLIA-waived urine drug screens, oral fluid drug and alcohol screens, and drug test controls.